



CERTIFICATE OF ANALYSIS

Product Name: CLC/Ipomoea CP Basal Medium
With Macronutrients, Micronutrients, & Vitamins
as described by Chee et al. (1992)

MADE IN THE USA

Catalog Number: CPP01

Storage: 2 to 8°C

Manufacture Date: APR-2014

Lot Number: 04140001

Expiration Date: 30-APR-2019

Test	Specifications	Results
Appearance	White to Yellow Powder	PASS
Solubility in Water	6.73 g/L in \leq 30 minutes	17 minutes
Moisture Content	\leq 2%	< 0.75%
Particle Size	\geq 80% through 100 mesh	\geq 80%
pH (<i>measured at 24°C</i>)	Record	3.93
Osmolality	Record mM/Kg	134 mM/Kg
Plant Tissue Culture of: <i>Hosta spp.</i> <i>Vitis spp.</i> <i>Solanum tuberosum</i> <i>Tradescantia zebrina</i>	Supports plant growth with no morphological abnormalities	PASS

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This Product is intended for Laboratory Use Only.

Prepared by:

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.

Kurt Rieben
Quality Control

28-JUN-2018