



CERTIFICATE OF ANALYSIS

Product Name: DCR Basal Salts Mixture

MADE IN THE USA

Catalog Number: DCP01

Storage: 2 to 8°C

Manufacture Date: MAR-2019

Lot Number: 03190010

Expiration Date: 31-MAR-2023

Test	Specifications	Results
Appearance	White to Yellow Powder	PASS
Solubility in Water	1.64 g/L in ≤ 30 minutes	14 minutes, with warming
Moisture Content	Record %	6.1%
Particle Size	Record % through 100 mesh	76%
pH (<i>measured at 42° C</i>)	Record	4.14
Osmolality	Record mM/Kg	24 mM/Kg
Plant Tissue Culture of: <i>Solanum tuberosum</i> <i>Hosta spp.</i> <i>Vitis spp.</i> <i>Tradescantia zebrina</i>	Supports plant growth with no morphological abnormalities	PASS

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This Product is intended for Laboratory Use Only.

Prepared by:

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.

Kurt Rieben
Quality Control

21-MAR-2019