



CERTIFICATE OF ANALYSIS

Product Name: Dulbecco's Modified Eagle's / Ham's F-12 Medium (1:1)
With L-Glutamine
Without Phenol Red

MADE IN THE USA

Catalog Number: DFL16

Storage: 2 to 8°C

Manufacture Date: MAR-2019

Lot Number: 03191021

Expiration Date: 30-SEP-2020

Test	Specifications	Results
Appearance	Clear, Light-Yellow Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH <i>(measured at ~22°C)</i>	7.1 - 7.4	7.36
Osmolality	Record mM/Kg	307 mM/Kg
Endotoxin	≤ 1.0 EU/mL	< 0.1 EU/mL
	<i>Endotoxin testing was performed using the LAL detection method.</i>	
Mycoplasma	Negative	NEGATIVE
	<i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>	
Animal Cell Culture		
Cell Growth Promotion	3 Passages	PASS
Viability	≥ 80%	> 90%
	<i>Medium was supplemented with 10% FBS, then tested using CHO cells grown in a monolayer in a 5% CO₂ atmosphere.</i>	

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.

Jacob Scovill
Quality Control

8-APR-2019