



CERTIFICATE OF ANALYSIS

Product Name: RPMI 1640 Medium
With L-Glutamine
Without L-Leucine

MADE IN THE USA

Catalog Number: RPL05

Storage: 2 to 8°C

Manufacture Date: FEB-2018

Lot Number: 02181007

Expiration Date: 31-AUG-2019

Test	Specifications	Results
Appearance	Clear, Red-Orange Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH (<i>measured at 24°C</i>)	6.90 - 7.40	7.26
Osmolality	264 - 294 mM/Kg	283 mM/Kg
Endotoxin	≤ 1.0 EU/mL	< 0.1 EU/mL
	<i>Endotoxin testing was performed using the LAL detection method.</i>	
Mycoplasma	Negative	NEGATIVE
	<i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay.</i>	
Animal Cell Culture		
Cell Growth Promotion	3 Passages	PASS
Viability	≥ 80%	> 90%
	<i>Medium was supplemented with 10% FBS and 50 mg/L L-Leucine, then tested using VERO cells grown in a monolayer in a 5% CO₂ atmosphere.</i>	

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for laboratory use only.

Prepared by:

Reviewed by:

Kurt Rieben
Quality Control

22-FEB-2018

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.